

Anti-LAMP2A Rabbit antibody



Catalog Number: 166766

Product name

Rabbit IgG

Anti-LAMP2A Rabbit antibody

Applications

Specificity

WB, IHC-P, FC, IP

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

LAMP2A Rabbit polyclonal antibody

IHC: 1/20-1/100

Preparation

FC: 1/200

Antigen: A synthetic peptide of human LAMP2A

IP: 1/20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage


Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Validations



Western blot detection of LAMP2A in 3T3 cell lysates using LAMP2A Rabbit pAb(1:1000 diluted). Predicted band size: 45kDa. Observed band size: 110kDa.