

Anti-CYP1B1 Rabbit antibody



Catalog Number: 166788

Product name

Anti-CYP1B1 Rabbit antibody

WB: 1/2000

IHC: 1/50

Specificity

Human, Mouse, Rat

FC: 1/50

Antibody description

CYP1B1 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human CYP1B1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, FC

Dilutions

Validations



Western blot detection of CYP1B1 in Rat Brain lysates using CYP1B1 Rabbit pAb(1:1000 diluted). Predicted band size:61kDa.Observed band size:61kDa.



Immunohistochemistry of CYP1B1 in paraffin-embedded Human breast cancer tissue using CYP1B1 Rabbit pAb at dilution 1/50