

# Anti-Phospho-EGFR (Tyr1092) Rabbit antibody



Catalog Number: 166791

## Product name

Anti-Phospho-EGFR (Tyr1092) Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse

## Applications

WB, IHC-P, ICC/IF

## Antibody description

Phospho-EGFR (Tyr1092) Rabbit polyclonal antibody

## Dilutions

WB: 1/1000

IHC: 1/50

## Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Tyr1092 of human EGFR

ICC/IF: 1/50

## Validations

## Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal



Western blot detection of EGFR (Phospho-Tyr1092) in C6,3T3,HeLa cell lysates using EGFR (Phospho-Tyr1092) Rabbit pAb(1:1000 diluted).Predicted band size:134kDa.Observed band size:175kDa.