

# Anti-MAP1LC3A Rabbit antibody



Catalog Number: 166794

## Product name

Anti-MAP1LC3A Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

MAP1LC3A Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human MAP1LC3A

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/2000


IHC: 1/20

ICC/IF: 1/20-1/100

FC: 1/20-1/100

IP: 1/20

## Validations



Western blot detection of MAP1LC3A in K562, Rat Brain, HeLa cell lysates using MAP1LC3A Rabbit pAb(1:1000 diluted). Predicted band size: 14kDa. Observed band size: 16kDa.