

Anti-ACE2 Rabbit antibody



Catalog Number: 166796

General Information

Product name:

Anti-ACE2 Rabbit antibody

Immunogen:

Applications:

WB, IHC-P, ICC/IF, IP

Specificity:

Human, Mouse, Rat

Antibody description:

ACE2 Rabbit polyclonal antibody

Properties

Preparation:

Antigen: A synthesized peptide derived from human ACE2

Formulation:

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage:

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality:

Polyclonal antibody

Ig Type:Rabbit IgG

Applications


WB: 1/1000-1/2000

IHC: 1/50-1/200


ICC/IF: 1/50-1/200

IP: 1/50

Validations



Immunohistochemical analysis of paraffin-embedded human kidney



Western blot analysis of ACE2 expression in Human kidney lysate.