

# Anti-FOXC1 Rabbit antibody



Catalog Number: 166799

## Product name

Anti-FOXC1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FOXC1 Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human FOXC1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/1000


IHC: 1/200

ICC/IF: 1/100

FC: 1/100

IP: 1/20

## Validations



Western blot detection of FOXC1 in 3T3 cell lysates using FOXC1 Rabbit pAb(1:1000 diluted). Predicted band size:57kDa. Observed band size:75kDa.