

Anti-RAC3 Rabbit antibody



Catalog Number: 166804

Product name

Anti-RAC3 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

RAC3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human RAC3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC-P

Dilutions

WB: 1/1000

IHC: 1/20

Validations



Western blot detection of RAC3 in K562,Rat Brain cell lysates using RAC3 Rabbit pAb(1:1000 diluted).Predicted band size:21kDa.Observed band size:21kDa.