

# Anti-DDX21 Rabbit antibody



Catalog Number: 166809

## Product name

Rabbit IgG

Anti-DDX21 Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF

Human

## Dilutions

## Antibody description

WB: 1/1000

DDX21 Rabbit polyclonal antibody

IHC: 1/20

## Preparation

ICC/IF: 1/100

Antigen: A synthetic peptide of human DDX21

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage


Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

## Validations



Western blot detection of DDX21 in HeLa cell lysates using DDX21 Rabbit pAb(1:1000 diluted). Predicted band size:87kDa.Observed band size:87kDa.