

Anti-SMARCC1/BAF155 Rabbit antibody



Catalog Number: 166827

Product name

Anti-SMARCC1/BAF155 Rabbit antibody

Specificity

Human, Rat

Antibody description

SMARCC1/BAF155 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human SMARCC1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, ICC/IF, FC, IP

Dilutions


WB: 1/1000-1/5000

ICC/IF: 1/50-1/200

FC: 1/20

IP: 1/20-1/50

Validations



Western blot detection of SMARCC1/BAF155 in K562,C6,Hela cell lysates using SMARCC1/BAF155 Rabbit pAb(1:1000 diluted).Predicted band size:123kDa.Observed band size:155kDa.