

Anti-HRH3/Histamine 3 Receptor Rabbit antibody



Catalog Number: 166831

Product name

Rabbit IgG

Anti-HRH3/Histamine 3 Receptor Rabbit antibody

Applications

Specificity

WB, ICC/IF, FC

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1/1000

HRH3/Histamine 3 Receptor Rabbit polyclonal antibody

ICC/IF: 1/50

Preparation


FC: 1/20-1/100

Antigen: A synthetic peptide of human HRH3/Histamine 3 Receptor

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of HRH3/Histamine 3 Receptor in K562, C6, 3T3 cell lysates using HRH3/Histamine 3 Receptor Rabbit pAb(1:1000 diluted). Predicted band size: 49kDa. Observed band size: 49kDa.

Clonality

Polyclonal

Ig Type