

Anti-PCK2 Rabbit antibody



Catalog Number: 166854

Product name

Anti-PCK2 Rabbit antibody

WB: 1/1000-1/5000

IHC: 1/100

Specificity

Human, Mouse, Rat

ICC/IF: 1/100

FC: 1/20

Antibody description

PCK2 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human PCK2

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC

Dilutions

Validations



Western blot detection of PCK2 in C6,3T3,HeLa cell lysates using PCK2 Rabbit pAb(1:1000 diluted).Predicted band size:71kDa.Observed band size:71kDa.



Immunohistochemistry of PCK2 in paraffin-embedded Human colon cancer tissue using PCK2 Rabbit pAb at dilution 1/50