

# Anti-MonoMethyl-Histone H3 (Lys9) Rabbit antibody



Catalog Number: 166857

## Product name

Anti-MonoMethyl-Histone H3 (Lys9) Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

MonoMethyl-Histone H3 (Lys9) Rabbit polyclonal antibody

## Preparation

Antigen: Synthetic Peptide of Histone H3 (Mono Methyl Lys9)

## Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/500-1000

## Validations



Western blot analysis of 1) HeLa, 2) Raw264.7, 3) Rat Testis Tissue with Histone H3 (Mono Methyl Lys9) Rabbit pAb diluted at 1:1,000.