

Anti-eIF1A Rabbit antibody



Catalog Number: 166872

General Information

Product name:

Anti-eIF1A Rabbit antibody

Immunogen:

Applications:

WB, ICC/IF, IP

Specificity:

Human, Mouse, Rat

Antibody description:

eIF1A Rabbit polyclonal antibody

Properties

Preparation:

Antigen: A synthetic peptide of human eIF1A

Formulation:

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage:

Store at -20°C. Stable for 12 months from date of receipt.

Clonality:

Polyclonal antibody

Ig Type:Rabbit IgG


Applications

WB: 1/2000-1/10000

ICC/IF: 1/20

IP: 1/20

Validations



Western blot detection of eIF1A in K562,C6,3T3,HeLa cell lysates using eIF1A Rabbit pAb(1:1000 diluted).Predicted band size:17kDa.Observed band size:17kDa.