

# Anti-Cleaved PARP1 Rabbit antibody



Catalog Number: 166876

## Product name

Anti-Cleaved PARP1 Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human

## Applications

WB

## Antibody description

Cleaved PARP1 Rabbit polyclonal antibody

## Dilutions

WB: 1/1000

## Preparation

Antigen: A synthetic peptide of human Cleaved PARP1

## Validations

## Formulation


Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal



Western blot detection of Cleaved PARP1 in K562 cell lysates using Cleaved PARP1 Rabbit pAb(1:1000 diluted). Predicted band size: 113kDa. Observed band size: 89kDa.