

# Anti-TMEM9B Rabbit antibody



Catalog Number: 166877

## Product name

Anti-TMEM9B Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

TMEM9B Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human TMEM9B

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000-1/2000

## Validations



Western blot detection of TMEM9B in 3T3 cell lysates using TMEM9B Rabbit pAb(1:1000 diluted).Predicted band size:23kDa.Observed band size:23kDa.