

Anti-GDF15 Rabbit antibody



Catalog Number: 166903

Product name

Anti-GDF15 Rabbit antibody

Specificity

Human

Antibody description

GDF15 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human GDF15

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, IP


Dilutions

WB: 1/1000

IHC: 1/20

IP: 1/20

Validations



Western blot detection of GDF15/Mic-1 in HeLa, MCF7 cell lysates using GDF15/Mic-1 Rabbit pAb(1:1000 diluted). Predicted band size: 34 kDa. Observed band size: 35, 13 kDa.