

Anti-Cleaved-Caspase-3 p17 Rabbit antibody



Catalog Number: 166905

Product name

Anti-Cleaved-Caspase-3 p17 Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Mouse, Human

Applications

WB

Antibody description

Cleaved-Caspase-3 p17 Rabbit polyclonal antibody

Dilutions

WB: 1/1000


Preparation

Antigen: A synthetic peptide of mouse Caspase-3

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of Caspase-3 in C6,3T3,HeLa cell lysates using Caspase-3 Rabbit pAb(1:1000 diluted).Predicted band size:32kDa.Observed band size:32/17kDa.

Clonality

Polyclonal