

Anti-PI3 Kinase p85 alpha Rabbit antibody



Catalog Number: 166906

Product name

Anti-PI3 Kinase p85 alpha Rabbit antibody

Specificity

Human

Antibody description

PI3 Kinase p85 alpha Rabbit polyclonal antibody

Preparation

Antigen: Recombinant Protein of PI 3-kinase p85 α

Formulation

PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1/1000-2000

Validations

94KD
66KD
45KD
35KD
26KD



Western blot analysis of HepG2 with PI3 Kinase p85u03b1 Rabbit pAb diluted at 1:2,000.