

Anti-BMAL1 Rabbit antibody



Catalog Number: 166907

Product name

Anti-BMAL1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

BMAL1 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human BMAL1

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/500-1/2000

Validations



Western blot analysis of BMAL1 expression in human fetal heart lysate.