

# Anti-N Cadherin Rabbit antibody



Catalog Number: 166908

## Product name

Anti-N Cadherin Rabbit antibody

## Specificity

Human

## Antibody description

N Cadherin Rabbit polyclonal antibody

## Preparation

Antigen: Recombinant protein of human N Cadherin

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M Nacl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of N Cadherin in HepG2 cell lysates using N Cadherin Rabbit pAb(1:1000 diluted).Predicted band size:100KDa.Observed band size:140KDa.