

# Anti-E Cadherin Rabbit antibody



Catalog Number: 166921

## Product name

Anti-E Cadherin Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Mouse, Rat

## Applications

WB, IP

## Antibody description

E Cadherin Rabbit polyclonal antibody

## Dilutions


WB: 1/1000

IP: 1/20

## Preparation

Antigen: Recombinant protein of mouse E Cadherin

## Validations



## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

Western blot detection of E Cadherin in Mouse brain lysates using E Cadherin Rabbit pAb(1:1000 diluted).Predicted band size:98(cleavages),135KDa.Observed band size: 80-120(cleavages),135KDa.

## Clonality

Polyclonal