

Anti-XIAP Rabbit antibody

Catalog Number: 166922

Product name

Anti-XIAP Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

XIAP Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human XIAP

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of XIAP in K562, Rat Brain, 3T3 cell lysates using XIAP Rabbit pAb(1:1000 diluted). Predicted band size: 57kDa. Observed band size: 53kDa.