

# Anti-LYVE1 Rabbit antibody



Catalog Number: 166924

---

## Product name

Anti-LYVE1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

LYVE1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human LYVE1

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/500-1/2000

## Validations



Western blot analysis of LYVE1 expression in A549 cell lysate.