

# Anti-APPBP1 Rabbit antibody



Catalog Number: 166926

## Product name

Anti-APPBP1 Rabbit antibody

## Specificity

Human, Rat

## Antibody description

APPBP1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human APPBP1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, FC, IP


## Dilutions

WB: 1/1000-1/2000

FC: 1/20

IP: 1/20

## Validations



Western blot detection of APPBP1 in SKOV-3 cell lysates using APPBP1 Rabbit pAb(1:1000 diluted). Predicted band size:60KDa. Observed band size:60KDa.