

# Anti-GSTA3 Rabbit antibody



Catalog Number: 166929

## Product name

Anti-GSTA3 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

GSTA3 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human GSTA3

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of GSTA3 in C6,3T3,HeLa cell lysates using GSTA3 Rabbit pAb(1:1000 diluted).Predicted band size:25kDa.Observed band size:25kDa.