

Anti-ATG12 Rabbit antibody



Catalog Number: 166941

Product name

Anti-ATG12 Rabbit antibody

Specificity

Human

Antibody description

ATG12 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Apg12

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

WB: 1/1000-1/5000

ICC/IF: 1/50-1/200

Validations



Western blot detection of Apg12 in K562,3T3,HeLa cell lysates using Apg12 Rabbit pAb(1:1000 diluted).Predicted band size:15kDa.Observed band size:55kDa.