

Anti-Annexin VI Rabbit antibody



Catalog Number: 166955

Product name

Anti-Annexin VI Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Annexin VI Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human Annexin VI

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

ICC/IF: 1/20

Validations



Immunocytochemistry of Annexin VI(green) in Hela cells using Annexin VI Rabbit pAb at dilution



1/50, and DAPI(blue)

Western blot detection of Annexin VI in K562,Rat Brain,C6,3T3 cell lysates using Annexin VI Rabbit pAb(1:1000 diluted).Predicted band size:76kDa.Observed band size:70kDa.