

# Anti-JAK2 Rabbit antibody



Catalog Number: 166957

## Product name

Anti-JAK2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

JAK2 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human JAK2

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, ICC/IF

## Dilutions


WB: 1/500-1/2000

ICC/IF: 1/50-1/200

## Validations



Western blot analysis of JAK2 expression in Jurkat



cell lysate.

Immunofluorescent analysis of Jurkat cells, using JAK2 Antibody .