

# Anti-Rel B Rabbit antibody



Catalog Number: 166960

## Product name

Anti-Rel B Rabbit antibody

WB: 1/5000

ICC/IF: 1/20

## Specificity

Human

FC: 1/20

IP: 1/20

## Antibody description

Rel B Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human Rel B

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, ICC/IF, FC, IP

## Dilutions

## Validations



Western blot detection of Rel B in HeLa cell lysates using Rel B Rabbit pAb(1:500 diluted). Predicted band size:62kDa.Observed band size:62kDa.



Western blot detection of Rel B in Raji cell lysates using Rel B Rabbit pAb(1:1000 diluted). Predicted band size:62KDa.Observed band size:62KDa.