

Anti-Phospho-YAP1 (Ser127) Rabbit antibody



Catalog Number: 166964

Product name

WB, IHC-P

Anti-Phospho-YAP1 (Ser127) Rabbit antibody

Dilutions

Specificity

WB: 1/5000-1/20000

Human, Mouse, Rat

IHC: 1/50-1/200

Antibody description

Phospho-YAP1 (Ser127) Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human YAP1

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality


Polyclonal

Ig Type


Rabbit IgG

Applications

Validations



Western blot analysis of Phospho-YAP1 (S127) expression in (1) HeLa cell lysate; (2) HeLa cell lysate treated with FBS+calyculin A.



Immunohistochemical analysis of paraffin-embedded human stomach, using Phospho-YAP1 (S127) Antibody.