

Anti-CDK5RAP3 Rabbit antibody



Catalog Number: 166976

Product name

WB: 1/1000

Anti-CDK5RAP3 Rabbit antibody

IHC: 1/20

Specificity

FC: 1/20

Human, Mouse

IP: 1/20

Antibody description

CDK5RAP3 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human CDK5RAP3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, FC, IP

Dilutions

Validations



Immunohistochemistry of CDK5RAP3 in paraffin-embedded Human breast cancer tissue using CDK5RAP3 Rabbit pAb at dilution 1/100



Western blot detection of CDK5RAP3 in K562, Hela cell lysates using CDK5RAP3 Rabbit pAb (1:1000 diluted). Predicted band size: 57 kDa. Observed band size: 70 kDa.