

Anti-Dlx4/BP1 Rabbit antibody



Catalog Number: 166986

Product name

Anti-Dlx4/BP1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Dlx4/BP1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Dlx4

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP

Dilutions


WB: 1/1000

IHC: 1/20-1/100

ICC/IF: 1/20-1/100

IP: 1/20

Validations



Western blot detection of Dlx4 in K562, Rat Brain, HeLa cell lysates using Dlx4 Rabbit pAb(1:1000 diluted). Predicted band size: 26kDa. Observed band size: 35kDa.