

Anti-hnRNP M1-M4 Rabbit antibody



Catalog Number: 166993

Product name

Anti-hnRNP M1-M4 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

hnRNP M1-M4 Rabbit polyclonal antibody

Preparation

Antigen: CPAQNEKRKEKNIKR

Formulation

10mM PBS with 0.02% sodium azide 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/500

Validations



Western blot analysis of extracts from HeLa cell lysates using hnRNP M Rabbit pAb (1:500 diluted). Predicted band size: 70kDa. Observed band size: 70kDa.