

Anti-CD8 alpha Rabbit antibody



Catalog Number: 167014

Product name

Anti-CD8 alpha Rabbit antibody

Specificity

Human

Antibody description

CD8 alpha Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human CD8

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/500-1/2000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200


Validations



Immunohistochemical analysis of paraffin-embedded human colon, using CD8 Antibody.



Western blot analysis of extracts from Raji cells using CD8 Rabbit pAb at 1:1000 dilution. Predicted band size: 32kDa. Observed band size:




Immunofluorescent analysis of Jurkat cells, using

Anti-CD8 alpha Rabbit antibody



Catalog Number: 167014

Western blot analysis of CD8 expression in Jurkat cell lysate.



CD8 Antibody .