

# Anti-TTF1 Rabbit antibody



Catalog Number: 167026

## Product name

Anti-TTF1 Rabbit antibody


ICC/IF: 1/50-1/200

FC: 1/50

## Validations

## Specificity

Human, Mouse, Rat



## Antibody description

TTF1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human TTF1


## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Immunohistochemical analysis of paraffin-embedded human thyroid, using TTF1 Antibody.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.



## Clonality

Polyclonal

Western blot analysis of TTF1 expression in 293 cell lysate.

## Ig Type

Rabbit IgG


## Applications

WB, IHC-P, ICC/IF, FC

## Dilutions

WB: 1/500-1/2000

IHC: 1/50-1/200



Immunofluorescent analysis of HeLa cells, using TTF1 Antibody .