

Anti-NLRP3 Rabbit antibody



Catalog Number: 167029

Product name

Anti-NLRP3 Rabbit antibody

Specificity

Mouse, Rat

Antibody description

NLRP3 Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human NALP3

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP, FC

Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200


ICC/IF: 1/50-1/200

IP: 1/20

Validations




Immunohistochemical analysis of paraffin-embedded human lung carcinoma, using NALP3



Antibody.

Immunofluorescent analysis of Raw264.7 cells,



using NALP3 Antibody.

Western blot analysis of NLRP3 expression in SH-SY5Y cell lysate.