

Anti-CDKN2A/p16INK4a Rabbit antibody



Catalog Number: 167036

Product name

Anti-CDKN2A/p16INK4a Rabbit antibody

Specificity

Human

Antibody description

CDKN2A/p16INK4a Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human p16 INK

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, IP, FC

Dilutions

WB: 1/500-1/2000


IHC: 1/50-1/200

ICC/IF: 1/50-1/200

IP: 1/30

FC: 1/100

Validations



Western blot detection of CDKN2A/p16INK4a in 293T, HeLa cell lysates using CDKN2A/p16INK4a Rabbit pAb(1:1000 diluted). Predicted band size: 17kDa. Observed band size: 17kDa.