

Anti-Carbonic Anhydrase 1/CA1 Rabbit antibody



Catalog Number: 167038

Product name

Anti-Carbonic Anhydrase 1/CA1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Carbonic Anhydrase 1/CA1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Carbonic Anhydrase I

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF


Dilutions

WB: 1/2000-1/10000


IHC: 1/50

ICC/IF: 1/200-1/1000

Validations



Immunohistochemistry of Carbonic Anhydrase I in paraffin-embedded Human colon cancer tissue using Carbonic Anhydrase I Rabbit pAb at dilution



Western blot detection of Carbonic Anhydrase 1/CA1 in Rat Brain lysates using Carbonic Anhydrase 1/CA1 Rabbit pAb(1:1000 diluted). Predicted band size:29kDa.Observed band size:29kDa.