

Anti-IL-1RA Rabbit antibody



Catalog Number: 167047

Product name

Anti-IL-1RA Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

IL-1RA Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human IL1RA

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/5000


IHC: 1/50

ICC/IF: 1/100

FC: 1/20

IP: 1/20

Validations



Western blot detection of IL-1RA in K562,3T3,HeLa cell lysates using IL-1RA Rabbit pAb(1:1000 diluted).Predicted band size:20kDa.Observed band size:20kDa.