

# Anti-FBP1 Rabbit antibody



Catalog Number: 167050

## Product name

Anti-FBP1 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

FBP1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human FBP1

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions


WB: 1/1000-1/5000

IHC: 1/20-1/50

## Validations



Immunohistochemistry of FBP1 in paraffin-embedded Human lung cancer tissue using FBP1



Rabbit pAb at dilution 1/20

Western blot detection of FBP1 in Rat Brain lysates using FBP1 Rabbit pAb(1:1000 diluted). Predicted band size: 37kDa. Observed band size: 37kDa.