

# Anti-MGST3 Rabbit antibody



Catalog Number: 167053

## Product name

Anti-MGST3 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

MGST3 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human MGST3

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, IP


## Dilutions

WB: 1/1000

IHC: 1/100

IP: 1/20

## Validations



Western blot detection of MGST3 in K562, HeLa cell lysates using MGST3 Rabbit pAb(1:1000 diluted). Predicted band size: 17kDa. Observed band size: 17kDa.