

Anti-NUMB Rabbit antibody



Catalog Number: 167063

Product name

Anti-NUMB Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

NUMB Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human NUMB

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, FC, IP

Dilutions


WB: 1/1000

IHC: 1/200

FC: 1/100

IP: 1/20

Validations



Western blot detection of NUMB in 3T3,C6,Rat Brain,MCF7 cell lysates using NUMB Rabbit pAb(1:1000 diluted).Predicted band size:71kDa.Observed band size:71kDa.