

Anti-TTF1 Rabbit antibody



Catalog Number: 167072

Product name

Anti-TTF1 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TTF1 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human TTF1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP

Dilutions

WB: 1/1000


IHC: 1/50

ICC/IF: 1/50

FC: 1/20-1/100

IP: 1/20

Validations



Western blot detection of TTF1 in 3T3 cell lysates using TTF1 Rabbit pAb(1:1000 diluted).Predicted band size:39kDa.Observed band size:39kDa.