

# Anti-Histone H2A.Z Rabbit antibody



Catalog Number: 167076

## Product name

Anti-Histone H2A.Z Rabbit antibody

## Ig Type

Rabbit IgG

## Specificity

Human, Mouse, Rat

## Applications

WB

## Antibody description

Histone H2A.Z Rabbit polyclonal antibody


## Dilutions

WB: 1/1000

## Preparation

Antigen: A synthetic peptide of human Histone H2A.Z

## Validations



## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

Western blot detection of Histone H2A.Z in K562,C6,3T3 cell lysates using Histone H2A.Z Rabbit pAb(1:1000 diluted).Predicted band size:14kDa.Observed band size:14kDa.