

Anti-NDUFAF5 Rabbit antibody



Catalog Number: 167082

Product name

Anti-NDUFAF5 Rabbit antibody

WB: 1/1000

IHC: 1/200

Specificity

Human, Mouse, Rat

ICC/IF: 1/100

FC: 1/50

IP: 1/20

Antibody description

NDUFAF5 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human NDUFAF5

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, FC, IP

Dilutions

Validations



Immunohistochemistry of NDUFAF5 in paraffin-embedded Human breast cancer tissue using



NDUFAF5 Rabbit pAb at dilution 1/50

Western blot detection of NDUFAF5 in 3T3 cell lysates using NDUFAF5 Rabbit pAb(1:1000 diluted). Predicted band size:39kDa.Observed band size:39kDa.