

# Anti-LAT Rabbit antibody



Catalog Number: 167086

---

## Product name

Anti-LAT Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

LAT Rabbit polyclonal antibody

## Preparation

Antigen: Synthesized peptide derived from LAT .  
at AA range: 190-270

## Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA  
and 0.02% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, ELISA

## Dilutions

WB: 1/500 - 1/2000

ELISA: 1/10000

## Validations



Western blot analysis of 293 using LAT antibody.  
Secondary antibody was diluted at 1:20000