

Anti-Phospho-AMPK alpha 2 (Thr172) Rabbit antibody



Catalog Number: 167089

Product name

Anti-Phospho-AMPK alpha 2 (Thr172) Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Phospho-AMPK alpha 2 (Thr172) Rabbit polyclonal antibody

Preparation

Antigen: A synthesized peptide derived from human AMPK alpha

Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/500-1/2000

Validations



Western blot analysis of Phospho-AMPK alpha (T172) expression in Mouse heart lysate.