

# Anti-Phospho-AMPK alpha 2 (Thr172) Rabbit antibody



Catalog Number: 167089

## Product name

Anti-Phospho-AMPK alpha 2 (Thr172) Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Phospho-AMPK alpha 2 (Thr172) Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human AMPK alpha

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/500-1/2000

## Validations



Western blot analysis of Phospho-AMPK alpha (T172) expression in Mouse heart lysate.