

# Anti-Rb Rabbit antibody



Catalog Number: 167104

## Product name

Rabbit IgG

Anti-Rb Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, ELISA

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/500 - 1/2000

Rb Rabbit polyclonal antibody

IHC: 1/100 - 1/300

## Preparation

IF: 1/200 - 1/1000

Antigen: The antiserum was produced against synthesized peptide derived from human Retinoblastoma. AA range:751-800

ELISA: 1/20000

## Validations

## Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type



Western blot detection of Rb in COS7 and C6 cell lysates using Rb Rabbit pAb (1:2000 diluted). Predicted band size:106KDa.Observed band size:110KDa.