

# Anti-GTPase HRAS Rabbit antibody



Catalog Number: 167106

## Product name

IgG

Anti-GTPase HRAS Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ELISA

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/500 - 1/2000

GTPase HRAS Rabbit polyclonal antibody

IHC: 1/100 - 1/300

## Preparation

ELISA: 1/20000

Antigen: The antiserum was produced against synthesized peptide derived from human RASH/RASK. AA range:1-50

## Validations

## Formulation

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.


## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type



Western blot detection of N/H/K-Ras in T47D,3T3,COS7 and HeLa cell lysates using N/H/K-Ras rabbit pAb (1:1000 diluted).Predicted band size:17kDa.Observed band size:17kDa.